HIV 1/2 STAT-PAK[®] DIPSTICK

A Qualitative Screening Test Kit for the Detection of Antibodies to HIV 1/2 in Human Serum, Plasma or Whole Blood

**FOR EXPORT ONLY
FOR *IN VITRO* DIAGNOSTIC USE
FOR PROFESSIONAL USE ONLY**

READ INSTRUCTIONS FOR USE CAREFULLY BEFORE PERFORMING TEST

INTENDED USE

The Chembio HIV 1/2 STAT-PAK DIPSTICK assay is a single use, immunochromatographic screening test which uses a cocktail of antigens to detect antibodies to HIV 1 and 2 in serum, plasma or whole blood. Reactive results are supportive evidence of exposure to HIV 1/2 and can be used to support a clinical diagnosis of HIV 1 or HIV 2. Nonreactive results, however, should not be used to exclude infection with HIV 1 or 2.

SUMMARY AND EXPLANATION

Discovered in 1983, Human Immunodeficiency Virus (HIV) is a retrovirus identified as the etiologic agent for Acquired Immunodeficiency Syndrome (AIDS) [1]. AIDS is characterized by changes in the population of T-cell lymphocytes that play a key role in the immune defense system. In the infected individual the virus causes a depletion of subpopulation of T-cells, called T-helper cells, which leaves these patients susceptible to opportunistic infections and certain malignancies. The major routes of transmission are sexual contact, contamination by blood or blood products and mother-to-newborn transmission [2-4].

Although there have been decreases in the rate of infection in certain countries, the number of persons infected with HIV globally has continued to increase. By the end of 2005 there were approximately 40.3 million people living with HIV/AIDS, up from about 37.5 million in 2003. Nearly 5 million people were newly infected with HIV/AIDS in 2005. In the same year more than 3 million died of AIDS-related illnesses; more than 500,000 of these were children [5].

The HIV virus consists of a genomic RNA molecule protected by a capsid and an envelope. The HIV envelope is the major target for humoral antibody response. The presence of the virus in patients causes the immune system to elicit the production of anti-HIV antibodies. The detection of these antibodies can be used as a diagnostic tool.

ELISA, Western Blots, PCR-based assays and various other test systems are currently available for HIV1/2 detection [6-10]. The Chembio HIV 1/2 STAT-PAK DIPSTICK assay is a rapid immunochromatographic test, which is simple and easy to use. The Chembio HIV 1/2 STAT-PAK DIPSTICK assay utilizes immobilized antigens for the detection of antibodies to HIV 1 and 2 in serum, plasma or whole blood with high degree of sensitivity and specificity.

PRINCIPLE OF TEST

The Chembio HIV1/2 STAT-PAK DIPSTICK assay employs a combination of antibody binding protein, which is conjugated to colloidal gold dye particles, and antigens to HIV1/2, which are bound to the membrane solid phase. The sample being tested and Running Buffer are applied to the sample pad. The Running Buffer facilitates the lateral flow of the specimen through the membrane and promotes the binding of antibodies to the antigens. If present, the antibodies bind to the gold conjugated antibody binding protein. In a reactive sample, the dye conjugated-immune complex migrates on the nitrocellulose membrane and is captured by the antigens immobilized in the TEST area producing a pink/purple line. In the absence of HIV 1/2 antibodies, there is no pink/purple line in the TEST area. The sample continues to migrate along the membrane and produces a pink/purple line in the CONTROL area containing immunoglobulin G antigens. This procedural control serves to demonstrate that specimens and reagents have been applied properly and have migrated through the device.

MATERIALS AVAILABLE AS AN ACCESSORY TO THE KIT

Chembio Rapid HIV 1/2 Antibody Test Kit Controls (Catalogue HIV104). Each package contains:

- HIV 1 Reactive Control
- HIV 2 Reactive Control
- Non Reactive Control
- Product Insert for the HIV104 Kit Controls

MATERIALS PROVIDED

Each kit contains the items to perform 30 tests:

- 30 HIV 1/2 STAT-PAK DIPSTICK test strips
- 1 Running Buffer bottle (8 ml)
- 30 Disposable Sample Loops (5 µL)
- 30 Backing Cards
- 1 Product Insert

MATERIALS REQUIRED BUT NOT PROVIDED

- Clock, watch or other timing device
- Automatic pipettor capable of delivering a 5 µL sample may be used in lieu of the disposable 5 µL sample loop supplied with the kit, for other than fingerstick specimens
- Disposable gloves
- Sterile Gauze (for fingerstick whole blood specimens)
- Sterile Alcohol swab (for fingerstick samples only)
- Biohazard disposal container
- Collection devices for samples other than fingerstick

WARNINGS

For *in-vitro* diagnostic use

1. Read the Product Insert completely before using this assay. Follow the instructions carefully as not doing so may result in inaccurate test results.
2. Use of the test kit with specimen types other than those specifically approved for use with this device may result in inaccurate test results.
3. This test should be performed at 18-30°C (64-86°F). If stored refrigerated, ensure that the test device is brought to operating temperature before performing testing.
4. If the test kit is stored at temperatures outside the storage temperature 8-30°C (46-86°F), or used outside the operating temperature 18-30°C (64-86°F), use the Kit Controls (Catalogue HIV104) to ensure proper performance of the test.
5. Individual(s) infected with HIV-1 and/or HIV-2 who is receiving highly active antiretroviral therapy (HAART) may produce negative results.
6. Sample loops are single use only. Discard sample loops following initial use.

STORAGE AND STABILITY

The HIV1/2 STAT-PAK DIPSTICK test strips should be stored at temperatures between 8-30°C in the original vial with desiccant. The Running Buffer should be stored at 8-30°C.

The vial containing the test strips should only be opened to remove dipsticks for testing and should be closed immediately and tightly thereafter.

NOTE: Do not use test kits beyond the expiration date imprinted on the box label and/or vial.

CAUTION: DO NOT FREEZE TEST KITS

PRECAUTIONS

Safety Precautions

1. Do not eat, drink or smoke in the area where specimens and kit reagents are handled. Avoid any contact between hands, eyes or mouth during specimen collection and testing.
2. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when handling patient specimens.
3. Dispose of all specimens and materials used in the test procedure in a biohazard waste container. Lancets should be placed in a puncture-resistant container prior to disposal. The recommended method of disposal of biohazard waste is autoclaving for a minimum of 1 hour at 121°C. Disposable materials may be incinerated. Liquid wastes may be mixed with appropriate chemical disinfectants. A freshly prepared solution of 10% bleach (0.5% solution of sodium hypochlorite) is recommended. Allow 60 minutes for effective decontamination.

NOTE: Do not autoclave solutions that contain bleach.

4. For additional information refer to: World Health Organization (WHO) Rapid HIV Tests: Guidelines for use in HIV testing and counseling services in resource-constrained Settings (2004) [11].

Handling Precautions

1. Do not use any test strips if the tube has been deformed.
2. Each device is for single use only.
3. Each sample loop is single use only and must be discarded following initial use.

4. Always check expiration date prior to testing.
5. Do not use the test beyond the expiration date printed on the tube containing the test strips.
6. Strips should be used immediately following removal from the storage tube.
7. Do not mix reagents from different lot numbers of kits.
8. Adequate lighting is required to read the test results.

SPECIMEN COLLECTION

The Chembio HIV 1/2 STAT-PAK DIPSTICK assay can be performed on whole blood, serum or plasma.

Fingerstick Whole Blood:

Prepare to perform the fingerstick blood collection procedure. Clean the finger of the person being tested with an antiseptic wipe. Allow the finger to dry thoroughly or wipe dry with a sterile gauze pad.

Using a sterile lancet, following laboratory procedure, prick the finger just off the center of the finger and wipe away the first drop with sterile gauze. Avoid squeezing the fingertip to accelerate bleeding as this may dilute the blood with excess tissue fluid.

Collect the sample from the second drop by touching the disposable sample loop provided to the drop of blood until the sample loop is full. Test immediately, following test procedure instructions.

Venous Whole Blood:

Draw blood following laboratory procedure for obtaining venous blood. Depending on use, collect sample in a tube containing citrate, heparin or EDTA. Be sure the tube of blood is well mixed before sampling. Dip the sample loop into the blood and allow it to fill. Follow test procedure instructions.

If tested the same day, venous whole blood may be kept at room temperature. Venous whole blood may be stored for up to 3 days between 2-8°C (36-46°F) before testing.

DO NOT FREEZE WHOLE BLOOD! Allow for refrigerated sample to reach room temperature and mix gently before testing.

Serum or Plasma:

Draw blood following laboratory procedures for obtaining serum or plasma specimens. Collect specimen in a clean tube not containing any anticoagulant (serum) or in a tube containing citrate, heparin or EDTA (plasma), following standard laboratory procedures for handling and collecting specimens. All specimens should be centrifuged and handled following local clinical or laboratory procedures. Remove the serum from the clot as soon as possible to avoid hemolysis.

Patient samples perform best when tested immediately after collection. If not to be tested immediately, specimens should be refrigerated following collection at 2-8°C and used within 3 days of collection. If testing within 3 days is not possible, the specimens should be frozen (-20°C or colder).

NOTE: If specimens are to be shipped, they should be packed in compliance with regulations covering the transportation of etiologic agents.


CONFIRMATION OF TEST PERFORMANCE

When the test is complete, you will observe a pink/purple line in the CONTROL area of the test strip on nonreactive as well as reactive samples. This serves as an internal control and gives confirmation of proper test performance. Pink/purple lines in both the TEST and CONTROL areas indicate reactive samples.


TEST PROCEDURE

If specimen to be tested is refrigerated, remove it from the refrigerator and allow it to come to room temperature prior to testing.


1. Remove the backing cards and tear off one segment for each test to be run. Enter patient or specimen ID in the space marked ID on the card.

2. Remove the HIV 1/2 STAT-PAK DIPSTICK test strip from the vial. The vial with the test strips should only be opened to remove dipsticks for testing and should be closed tightly immediately thereafter.
3. Being careful not to touch the membrane, turn the test strip so the green tape is facing away and remove the red liner from the adhesive strip(s) on the back of the test strip. Discard the red liner.

4. Place the test strip onto the backing card in the marked space with the green tape facing up and the arrows on the tape facing in the same direction as the arrows on the backing card. Place the backing card and dipstick assembly on a clean, flat surface.

5. SPECIMEN COLLECTION

A. Fingertick Whole Blood

- i. Following laboratory procedure, prick the finger with a lancet and wipe away the first drop of blood.
- ii. Collect the sample from the second drop by touching the 5µL sample loop provided to the drop of blood. (Proceed to Step 6)


B. Serum, Plasma, Venous Whole Blood

Touch the 5µL sample loop provided to the sample to fill the circular opening of the loop with sample.


6. Holding the sample loop vertically, touch it to the sample pad in the center of the SAMPLE (S) area of the dipstick to dispense ~5 µL of sample (serum, plasma or whole blood) onto the sample pad.


- Invert the Running Buffer bottle and hold it vertically (not at an angle) over the sample area. Add 3 drops (~ 105 µL) of buffer slowly, dropwise, onto the SAMPLE (S) area being sure each drop is absorbed by the sample pad before adding the next.

NOTE: In the case that liquid does not start to flow onto the membrane within 2 minutes, add one additional drop of Running Buffer.

Control Line
Test Line

- Read the test result between 15 and 20 minutes after the addition of the Running Buffer. In some cases a test line may appear in less than 15 minutes however, 15 minutes are needed to report a nonreactive result. Read results in a well-lit area.

Do not read results after 20 minutes.

NOTE: Discard the used sample loop, test strip, backing card and gloves into a biohazard waste container.

QUALITY CONTROL

A pink/purple line should always appear in the CONTROL area if the test has been performed correctly and the device is working properly. It serves as an internal procedural control.

Good Laboratory Practice (GLP) recommends the use of control materials along with the test samples to ensure proper performance of the test kit. Reactive and nonreactive serum or plasma based commercial controls should be used for this purpose. Use controls, as per the TEST PROCEDURE instructions of this insert. Please note that certain commercial controls designed for ELISA may not perform properly with the Chembio HIV 1/2 STAT-PAK DIPSTICK assay. Serum or plasma based controls for example from BBI (Seracare/Boston Biomedica Inc, USA) or NABI (Nabi Diagnostics, USA) can be used.

RUN CONTROLS UNDER THE FOLLOWING CIRCUMSTANCES:

- Each new operator prior to performing tests on patient specimens,
- When opening a new kit lot,
- Whenever a new shipment of Kits is received,
- If the temperature of the test storage area falls outside of 8-30°C (46-86°F),
- If the temperature of the testing area falls outside of 18-30°C (64-86°F),
- At periodic intervals as indicated by the user facility.

INTERPRETATION OF RESULTS


Reactive Result

Two pink/purple lines, one in the TEST area and one in the CONTROL area indicate a reactive result. The line in TEST area may look different from the line in the CONTROL area. The intensity of the line in the TEST area will vary with the concentration of specific antibodies, from barely visible to very dark.

NOTE: Even a very faint line in the TEST area should be considered reactive. It is recommended that a reactive result be confirmed by Western Blot or IFA in accordance with the World Health Organization (WHO) recommendations. This test is suitable for use in multi-test algorithms designed for the statistical validation of rapid HIV test results. When multiple rapid HIV tests are available, this test should be used in appropriate multi-test algorithms.

Nonreactive Result


One pink/purple colored line in the CONTROL area, with no colored line in the TEST area indicates a nonreactive result. A nonreactive result after 15 minutes indicates that there are no HIV antibodies in the patient sample, but this result does not exclude HIV infection.

*Red Cell Line occurs on rare occasions with blood specimens, but has no impact on test results. Such a line should not be interpreted as a positive test result

Invalid

A pink/purple line should always appear in the CONTROL area, no matter if the TEST LINE appears or not. If there is no distinct pink/purple line visible in the CONTROL area, the test is invalid. It is recommended that the test be repeated with a new device.

LIMITATIONS OF THE PROCEDURE

1. The HIV 1/2 STAT-PAK DIPSTICK procedure and the interpretation of the results must be followed closely. It is a screening test designed for detecting antibodies against HIV 1 and 2 in human serum, plasma or whole blood. Results from testing other body fluids or of pooled serum or plasma should not be used.
2. Venous whole blood must be collected with an anticoagulant depending on use. For serum samples collect blood without anticoagulant.
3. Do not open the tube containing test strips until just prior to use. Close it immediately after removing required number of test strips.
4. Do not use kit contents beyond their labeled expiration date.
5. Ensure finger is completely dry before performing fingerstick.
6. Read results in a well lit area.
7. Reactive results must be confirmed with Western Blot or IFA according to WHO recommendations. The clinical evaluation of the patient's situation should be performed before a final diagnosis is made. A reactive result obtained with a single rapid test should not be used to diagnose HIV/AIDS infection even if antibodies to HIV 1 and 2 are present. A nonreactive result at any time does not preclude the possibility of infection with HIV 1/2.
8. An individual infected with HIV-1 and/or HIV-2 who is receiving highly active antiretroviral therapy (HAART) may produce false negative results.
9. This assay has not been evaluated for newborn screening, cord blood specimens, or for individuals less than 13 years of age.

EXPECTED RESULTS

This is a qualitative test for the detection of antibodies to HIV 1 and 2 in whole blood, serum or plasma. As described in the PERFORMANCE CHARACTERISTICS section below, the sensitivity of the Chembio HIV 1/2 STAT-PAK DIPSTICK assay was found to be substantially equivalent to EIA and Western Blot tests when tested on selected performance panels.

PERFORMANCE CHARACTERISTICS

In-house studies demonstrate that the sensitivity of the Chembio HIV 1/2 STAT-PAK DIPSTICK assay is substantially equivalent to EIA and Western Blot tests when tested on BBI performance panels PRZ204, Anti-HIV 1/2 Combo Performance Panel, and PRB204 Mixed Titer Performance Panel and WWRB302 and WWRB303 Worldwide Performance Panels.

Performance data were also obtained in three external studies using the alternate test procedure in a tube. In a limited evaluation of 770 samples by the World Health Organization (WHO), the final sensitivity was 99.0% and the final specificity was 100% compared to the reference assays. In this study, none of the results were recorded as indeterminate. Evaluations in Zambia with 236 specimens, (125 reactive / 111 non-reactive) and in Uganda with 503 specimens (90 reactive / 413 non-reactive) both resulted in 100% sensitivity and 100% specificity for the Chembio HIV 1/2 STAT-PAK DIPSTICK assay.

PRECISION

Intraassay

Within run precision was determined by using 10 replicates of two specimens containing different levels of HIV 1/2 antibodies. The non reactive and reactive results were correctly identified 100% of the time.

Interassay

Between run precision was determined by using the same two specimens in 10 different replicates from three different lots of test devices. Nonreactive and reactive results were observed 100% of the time.

Cross Reactivity and Interference

No cross reactivity was observed from Hepatitis B, Rheumatoid Factor (~80 IU/ml) and hCG (500 mIU/ml) with the HIV 1/2 STAT-PAK DIPSTICK assay. In addition, no interference from bilirubin, hemoglobin or triglycerides was observed.

REFERENCES

1. Essex, M. (1999) Human immunodeficiency viruses in the developing world. *Adv Virus Res* 53 : 71-88.
2. Kanki, P.J., Hopper, J.R. and Essex, M. (1987) The origins of HIV-1 and HTLV-4/HIV/2. *Ann N Y Acad Sci* 511 : 370-375.
3. Nicoll, A., Gill, O.N. (1999) The global impact of HIV infection and disease. *Commun Dis Publ Health* 2 : 85-95.
4. Valdiserri R.O., Holtgrave, D.R., West, G.R. (1999) Promoting early diagnosis and entry into care. *AIDS* 13 : 2317-2330.
5. UNAIDS/WHO Press Release (2005) HIV Infection Rates Decreasing in Several Countries but Global Number of People Living with HIV Continues to Rise. Available at www.unaids.org/epi/2005/doc/docs/PR_EpiUpdate_Nov05en.pdf : 1-3.
6. Essex, M., Kanki, P. J., Marlink, R., et al. (1990) Antigenic characterization of the human immunodeficiency viruses. *J Am Acad Dermatol* 22 : 1206-1210.
7. Essex, M., McLane, M.F., Lee, T.H. et al. (1983) Antibodies to cell membrane antigens associated with human T-cell leukemia virus in patients with AIDS. *Science* 220 : 859-862.
8. Gallo, R.C., Saluahuudin, S.Z., Popovic, M., et al. (1984) Frequent detection and isolation of cytopathic retroviruses (HTLV-III) from patients with AIDS and at risk for AIDS. *Science* 224 : 500-503.
9. Kenealy, W., Reed, D., Cybulsky, R., et al. (1987) Analysis of human serum antibodies to human immunodeficiency virus (HIV) using recombinant ENV and GAG antigens. *AIDS Res Human Retrovir* 3 : 95-105.
10. Kovacs, A., Xu, J., Rasheed, S., et al. (1995) Comparison of a rapid non-isotopic polymerase chain reaction assay with four commonly used methods for the early diagnosis of human Immunodeficiency Virus type 1 infection in neonates and children. *Pediatr Infect Dis J* 14 : 948-954.
11. World Health Organization (WHO): Rapid HIV Tests: Guidelines for Use in HIV Testing and Counseling Services in Resource – Constrained Settings. (2004).

SYMBOL LEGEND

	CONSULT THE MANUAL BEFORE USE
	CAUTION, CONSULT ACCOMPANYING DOCUMENTS.
	DO NOT REUSE
	FOR USE WITHIN TEMPERATURE LIMITS
	IN VITRO DIAGNOSTIC MEDICAL DEVICE
	BATCH CODE
	PRODUCT CATALOG NUMBER
	MANUFACTURERS IDENTIFICATION
	DATE OF MANUFACTURE
	USE BY DATE

FOR MORE INFORMATION, CONTACT

 CHEMBIO DIAGNOSTIC SYSTEMS, INC.
3661 HORSEBLOCK ROAD
MEDFORD, NY 11763 USA

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